

Download Biopure Case Analysis

Ideal for the quantitation of proteins, peptides and non-sequence specific total DNA. The production of biopharmaceuticals, such as monoclonal antibodies, recombinant proteins, vaccines, and nucleic acids requires special techniques and quality control methods. As you may know, there are several types of water purity available depending on what you're tackling at the bench, from Type III for general use such as rinsing out your beakers, all the way up to Type I+ for sensitive applications like the intriguingly sounding 'graphite furnace atomic absorption spectrometry' (GF-AAS). Mr. Leibensperger has 18 years of experience in Biopharmaceutical Manufacturing and Operations Management. He was an original team member with Biopure Corp. that helped bring Oxyglobin® and Hemopure® production to full scale and an integral part of achieving FDA and EMA GMP approval for manufacturing facilities. Enzo Biochem Inc. stock price, stock quotes and financial overviews from MarketWatch.